




CHUNLI ORTHOPEDICS GLOBAL INTELLIGENCE



# TITANIUM PLATE WITH LOOP

Beijing Chunlizhengda Medical Instruments Co., Ltd



**Beijing Chunlizhengda Medical Instruments Co., Ltd**

2007 Gaohejiefang Plaza, No.98 South of East 3rd Ring Road, Chaoyang Dist. Beijing (100021)

Mobile Hotline:

- +86 13810479615 (English)
- +86 18516853658 (Spanish)
- +86 13520826455 (Russian)



Company Wechat

[www.clzd.com](http://www.clzd.com)

# TITANIUM PLATE WITH LOOP


**SCOPE OF APPLICATION:** It is suitable for the fixation of ligament and muscle or bone during orthopedic reconstruction.

Specification	Description
TA/8	Oval plate (removable with openings at both ends)
TB/3.5,3、 4	Distal traction (adjustable)
TC/3.5,3、 4	Proximal traction (adjustable)
TD/8	Four-leaf clover shift (removable with openings at both ends)
TE/14	Round titanium plate (removable with openings at both ends)
TF/3、 4	Coil (adjustable only with coils)


TA



TB



TC



TD



TE



## PRODUCT FEATURES

- 1 Complete range, including a single plate with openings at both ends, a single coil with adjustable length, and a combination of plate and coil with adjustable coil length, are suitable for a variety of clinical demand.
- 2 The traction lines could be offered for the combined product to help the physician pull the graft indirectly for adjustment.
- 3 Products do not need to turn the loop line, reduce the use of steps.
- 4 The coil can pass through the bone canal of small diameter, the length is adjustable, can be suitable for different patients, and the coil end is suitable for various types of titanium plate.
- 5 Coils or a combination of products with coils allow maximum contact with bone.
- 6 Each titanium plate has grooves that allow stitches and knots to be held in place without causing noticeable bumps

TF

