


(BZ) NO NEEDLE SUTURE SPECIFICATIONS

Suture line no	Texture/Form	Color	Length specification
3-0#	Ultra high molecular weight polyethylene braided suture	White Blue Blue-white	45cm-130cm One gauge at 5cm intervals
2-0/T#			
2-0#			
0#			
1#			
2#			
3、4#			
5#			

CHUNLI ORTHOPEDICS GLOBAL INTELLIGENCE

CHUNLi 



CHUNLi 

Beijing Chunlizhengda Medical Instruments Co., Ltd

2007 Gaoheyanfeng Plaza, No.98 South of East 3rd Ring Road, Chaoyang Dist. Beijing (100021)

Mobile Hotline:

- +86 13810479615 (English)
- +86 18516853658 (Spanish)
- +86 13520826455 (Russian)



Company Wechat

www.clzd.com

NON-ABSORBABLE SURGICAL SUTURES

Beijing Chunlizhengda Medical Instruments Co., Ltd

NON-ABSORBABLE SURGICAL SUTURES

(DZ)SUTURE WITH NEEDLE SPECIFICATIONS

● Product composition

Non-absorbable surgical sutures are woven with ultra-high molecular weight polyethylene yarns conforming to YY/T 1431-2016, with some sutures with needles, and the stitches are made of 12Cr18Ni9 stainless steel material conforming to YY/T 0043-2016.

● Scope of application

Non-absorbable surgical suture products are suitable for suturing or ligation of soft tissues.

● Product features

- 1 Suture is woven with UHMWPE yarn with minimal tissue reaction and good biocompatibility;
- 2 Providing blue, white, blue-white color options for easy operation identification and suture management;
- 3 Providing a variety of specifications of different stitches and stitching needle combinations, complete product specifications, can meet the various needs of customers;
- 4 Suture surface is smooth, dry and uniform, no knot, colorfast, supple, wear-resistant, excellent knotting performance and tensile strength far beyond the standard requirements;
- 5 The suture needle is made of medical stainless steel, with excellent material, good elasticity and toughness, and sharp puncture. The needle and thread connection is smooth, the needle tail is smooth after connection, and the connection strength is high.



Suture line no	Specifications and models (Taking white,100cm, single needle as an example)	Suture shape (Single needle)	○Round needle △Cutting needle	Stitches arc Diameter× Chord length	Number of threads	Suture color	Suture strength									
3-0#	3-0#,100cm, White, ○3/8-6×13, S		○	3/8 6×13	(S) Single needle	White Blue Blue-	45cm- 130cm One size by 5cm interval									
	3-0#,100cm, White, ○3/8-6×14, S		○	3/8 6×14												
	3-0#,100cm, White, △1/2-7×17, S		△	1/2 7×17												
	3-0#,100cm, White, ○3/8-7×20, S		○	3/8 7×20												
	3-0#,100cm, White, ○3/8-7×17, S		○	3/8 7×17												
2-0#	2-0#,100cm, White, ○1/2-7×17, S		○	1/2 7×17				(D) Double needle	white							
	2-0#,100cm, White, ○3/8-7×17, S		○	3/8 7×17												
	2-0#,100cm, White, △3/8-8×20, S		△	3/8 8×20												
0#	0#,100cm, White, ○1/2-9×14, S		○	1/2 9×14							(D) Double needle	white				
	0#,100cm, White, ○1/2-9×24, S		○	1/2 9×24												
	0#,100cm, White, ○1/2-8×20, S		○	1/2 8×20												
1#	1#,100cm, White, △1/2-10×28, S		△	1/2 10×28										(D) Double needle	white	
	1#,100cm, White, ○1/2-10×34, S		○	1/2 10×34												
2#	2#,100cm, White, ○1/2-11×20, S		○	1/2 11×20	(D) Double needle	white										
	2#,100cm, White, ○1/2-11×24, S		○	1/2 11×24												
	2#,100cm, White, ○1/2-11×28, S		○	1/2 11×28												
	2#,100cm, White, ○1/2-11×34, S		○	1/2 11×34												
5#	5#,100cm, White, △1/2-12×30, S		△	1/2 12×30				(D) Double needle	white							
	5#,100cm, White, ○1/2-12×30, S		○	1/2 12×30												
	5#,100cm, White, ○1/2-12×34, S		○	1/2 12×34												
	5#,100cm, White, ○1/2-12×40, S		○	1/2 12×40												